

Package ‘dPCP’

June 23, 2026

Title Automated Analysis of Multiplex Digital PCR Data

Version 2.0.3

Description The automated clustering and quantification of the digital PCR data is based on the combination of 'DBSCAN' (Hahsler et al. (2019) <[doi:10.18637/jss.v091.i01](https://doi.org/10.18637/jss.v091.i01)>) and 'c-means' (Bezdek et al. (1981) <[doi:10.1007/978-1-4757-0450-1](https://doi.org/10.1007/978-1-4757-0450-1)>) algorithms.

The analysis is independent of multiplexing geometry, dPCR system, and input amount.

The details about input data and parameters are available in the vignette.

License MIT + file LICENSE

Encoding UTF-8

Depends R (>= 4.0.0)

Imports cluster, dbscan, e1071, exactci, ggplot2, ggpubr, graphics, raster, rlist, scales, shiny, shinyjs, stats, stringr, utils

Suggests knitr, rmarkdown, testthat

VignetteBuilder knitr

URL <https://github.com/alfodefalco/dPCP>

BugReports <https://github.com/alfodefalco/dPCP/issues>

Config/roxygen2/version 8.0.0

NeedsCompilation no

Author Alfonso De Falco [aut, cre] (ORCID: <<https://orcid.org/0000-0001-9985-9920>>),
Michel Mittelbronn [ctb],
Christophe Olinger [ctb],
Daniel Stieber [ctb],
Laboratoire national de santé [cph]

Maintainer Alfonso De Falco <alfonsodefalco90@gmail.com>

Repository CRAN

Date/Publication 2026-06-23 05:10:49 UTC

Contents

centers_data	2
cmeans_clus	3
dbscan_combination	5
dPCP	6
export_csv	9
manual_correction	10
rain_reclus	11
read_reference	12
read_sample	14
read_sampleTable	15
reference_dbscan	16
replicates_quant	18
report_dPCP	19
target_quant	20
Index	22

centers_data	<i>Prediction of clusters centroid position</i>
--------------	---

Description

This function calculates the coodintaes of all clusters centroid.

Usage

```
centers_data(sample.subquality, sample.table, referenceDB)

## S3 method for class 'centers_data'
plot(x, ..., sample = "all")
```

Arguments

- sample.subquality an object of class read_sample, inherited from [read_sample](#).
- sample.table object of class sample_table, inherited from [read_sampleTable](#).
- referenceDB an object of class reference_dbscan, inherited from [reference_dbscan](#)
- x an object of class centers_data
- ... Arguments to be passed to methods
- sample 'all' to show all samples, or a numeric vector indicating the row number of samples in the sample table.

Value

An object of class `centers_data` containing a sublist for each sample. Each sublist has the following components:

<code>quality</code>	quality threshold used in <code>read_sample</code> .
<code>reference</code>	reference ID.
<code>centers</code>	a data frame with the centroids coordinates.
<code>data</code>	a data frame with the fluorescence intensities.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
samp <- read_sample(sample.table, system = "bio-rad",
                   file.location = fileLoc)

#Reference DBSCAN clustering
dbref <- reference_dbscan(ref, sample.table, save.template = FALSE)

#Predict position of clusters centroid from reference DBSCAN results
cent <- centers_data(samp, sample.table, dbref)

plot(cent, sample = "all")

## End(Not run)
```

cmeans_clus

*Cluster analysis with fuzzy c-means algorithm***Description**

This function carries out the c-means cluster analysis, using the centroids position as initial values for cluster centers.

Usage

```
cmeans_clus(centers.data)

## S3 method for class 'cmeans_clus'
plot(x, ..., sample = "all", color.blind = FALSE)
```

Arguments

<code>centers.data</code>	an object of class <code>centers_data</code> , inherited from centers_data .
<code>x</code>	an object of class <code>cmeans_clus</code>
<code>...</code>	Arguments to be passed to methods
<code>sample</code>	'all' to show all samples, or a numeric vector indicating the row number of samples in the sample table.
<code>color.blind</code>	logical. If TRUE colors optimized for colorblind readers are used.

Value

An object of class `cmeans_clus` containing a sublist for each sample. Each sublist has the following components:

<code>quality</code>	quality threshold used in read_sample .
<code>reference</code>	reference ID.
<code>centers</code>	a data frame with the centroids coordinates.
<code>data</code>	a data frame with the fluorescence intensities and clusters name.
<code>membership</code>	a matrix with the membership values of the data elements to the clusters. See also cmeans

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
samp <- read_sample(sample.table, system = "bio-rad",
```

```

file.location = fileLoc)

#Reference DBSCAN clustering
dbref <- reference_dbscan(ref, sample.table, save.template = FALSE)

#Predict position of clusters centroid from reference DBSCAN results
cent <- centers_data(samp, sample.table,dbref)

#Fuzzy c-means clustering
cmclus <- cmeans_clus(cent)

plot(cmclus, sample = "all")

## End(Not run)

```

dbscan_combination *Test eps and minPts combinations for DBSCAN analysis*

Description

This function tests all combinations of eps and minPts for DBSCAN analysis of reference samples indicated in refID. The results are represented in scatterplots exported to a pdf file.

Usage

```

dbscan_combination(
  refID,
  system = NULL,
  file.location = ".",
  reference.quality = 0.5,
  eps = c(120, 150, 180, 200),
  minPts = c(20, 50, 80, 100)
)

```

Arguments

refID	a string or a character vector of chipID (Thermo Fisher) or the complete file name with the extension (Bio-Rad) of reference sample(s) to be analysed.
system	character. The name of digital PCR system used to generate the data. It must be either Thermo Fisher or Bio-Rad. Abbreviations are also accepted.
file.location	character. Full path name to reference and sample files location. The default corresponds to the working directory, (getwd). Tilde expansion (see (path.expand)) is performed.
reference.quality	numeric. Between 0 and 1. Quality threshold to subset the data (just for Thermo Fisher). If different thresholds have to be applied to various reference samples, a vector of the same length of refID has to be provided.

eps	a numeric vector of values to be tested. Maximum distance between elements within a cluster in a DBSCAN analysis. See also dbscan .
minPts	a numeric vector of values to be tested. Number of minimum elements to assemble a cluster in a DBSCAN analysis. See also dbscan .

Value

A pdf file containing the scatterplots of DBSCAN analysis performed with all combinations of eps and minPts. Each reference generates a different pdf file.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

dbscan_combination("dilution20200313_B01_Amplitude.csv",
                  file.location = fileLoc, system = "bio-rad",
                  eps = c(150, 160, 180, 190), minPts = c(80, 100, 120))

unlink("dilution20200313_B01_Amplitude.pdf")

## End(Not run)
```

dPCP

Automated analysis of digital PCR data

Description

This function carries out the automated clustering of digital PCR data.

Usage

```
dPCP(
  file,
  system = NULL,
  file.location = ".",
  reference.quality = 0.5,
  sample.quality = 0.5,
  eps = 200,
  minPts = 50,
  save.template = FALSE,
  rain = TRUE,
```

```

    QC.reference = FALSE,
    partition.volume = NULL
)

## S3 method for class 'dPCP'
plot(
  x,
  ...,
  sample = "all",
  reference = "all",
  type = "dPCP",
  color.blind = FALSE
)

```

Arguments

<code>file</code>	character. The name or the path of csv file to be read. If it does not contain an absolute path, the file name is relative to the current working directory, (getwd).
<code>system</code>	character. The name of digital PCR system used to generate the data. It must be either Thermo Fisher or Bio-Rad. Abbreviations are also accepted.
<code>file.location</code>	character. Full path name to reference and sample files location. The default corresponds to the working directory, (getwd). Tilde expansion (see (path.expand)) is performed.
<code>reference.quality</code>	numeric. Between 0 and 1. Quality threshold to subset the data. If different thresholds have to be applied to various reference samples, a vector of the same length of number of reference samples has to be provided. Used only when the system is Thermo Fisher.
<code>sample.quality</code>	numeric. Between 0 and 1. Quality threshold to subset data. If different thresholds have to be applied to various samples, a vector of the same length of number of samples has to be provided. Used only when the system is Thermo Fisher.
<code>eps</code>	numeric. Input parameter for the DBSCAN algorithm. It represents the maximum distance between the elements within a cluster. See also dbscan . If different values have to be applied to various reference samples, a vector of the same length of number of reference samples has to be provided.
<code>minPts</code>	numeric. Input parameter for the DBSCAN algorithm. It represents the number of minimum elements to assemble a cluster. See also dbscan . If different values have to be applied to various reference samples, a vector of the same length of number of reference samples has to be provided.
<code>save.template</code>	logical. If TRUE a template of DBSCAN analysis of reference samples is saved. When system is Thermo Fisher, <code>save.template</code> can be also a character vector indicating the chipID.
<code>rain</code>	logical. If TRUE the rain analysis is carried out.
<code>QC.reference</code>	logical. If TRUE the fraction of rain elements in the reference samples is calculated. Warning messages are displayed when the percentage of rain is high.

<code>partition.volume</code>	numeric. This parameters is taken into account when the parameter 'system' is set on Other. Indicate the partion volume in microliters spcific to the digital PCR system.
<code>x</code>	an object of class dPCP
<code>...</code>	Arguments to be passed to methods
<code>sample</code>	'all' to show all samples, or a numeric vector indicating the row number of samples in the sample table.
<code>reference</code>	'all' to show all reference samples, or a character vector with chip ID (Thermo Fisher) or the file name (Bio-rad) of reference samples to be showed.
<code>type</code>	string. Type of plot to be showed. Available plots: 'reference dbscan', 'centers', 'cmeans', 'rain', 'dPCP'. @param color.blind logical. If TRUE colors optimized for colorblind readers are used.
<code>color.blind</code>	logical. If TRUE colors optimized for colorblind readers are used.

Value

An object of class dPCP containing the following components:

<code>referenceDB</code>	an object of class reference_dbscan.
<code>samples</code>	a list of samples. Each sample sublist contains the information about the cluster analysis.
<code>results</code>	an object of class replicates_quant.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#dPCP analysis
results <- dPCP(sampleTable, system = "bio-rad", file.location = fileLoc,
               eps = 200, minPts = 50, save.template = FALSE, rain = TRUE,
               QC.reference = FALSE)

plot(results, sample = 1, type = "dPCP")

## End(Not run)
```

export_csv	<i>Export dPCP analysis results to a csv file</i>
------------	---

Description

This function exports dPCP analysis results to a csv file.

Usage

```
export_csv(data, filename)
```

Arguments

data	an object of class dPCP, target_quant or replicates_quant.
filename	character. File name (no extension) for csv and pdf files to create on disk.

Value

A csv file with the information and results of dPCP analysis.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata",package = "dPCP")

#dPCP analysis
results <- dPCP(sampleTable, system = "bio-rad", file.location = fileLoc,
               eps = 200, minPts = 50, save.template = FALSE,
               rain = TRUE)

export_csv(results, filename = "dPCRproject_1")

## End(Not run)
```

manual_correction	<i>Manual correction of dPCP cluster analysis</i>
-------------------	---

Description

This function builds an interactive app to manually correct the dPCP cluster analysis.

Usage

```
manual_correction(
  data,
  filename,
  save.plot = FALSE,
  format = "png",
  dpi = 300,
  color.blind = FALSE
)
```

Arguments

data	an object of class dPCP, inherited from dPCP .
filename	character. File name (no extension) for csv and pdf files to create on disk.
save.plot	logical. If TRUE the plots are exported to a file.
format	a string indicating the file format for the export. Available formats: 'eps', 'ps', 'tex', 'pdf', 'jpeg', 'tiff', 'png', 'bmp', 'svg', 'wmf'.
dpi	numeric. Image resolution.
color.blind	logical. If TRUE colors optimized for colorblind readers are used.

Value

A Shiny session.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#dPCP analysis
results <- dPCP(sampleTable, system = "bio-rad", file.location = fileLoc,
               eps = 200, minPts = 50, save.template = FALSE,
               rain = TRUE)
```

```

manual_correction(results, filename = "manual_dPCR", save.plot = FALSE)

## End(Not run)

```

rain_reclus

Identification and clustering of "rain" data

Description

This function identifies the "rain" elements and re-clusters them using the Mahalanobis distance. Each "rain" element is assigned to the cluster whose Mahalanobis distance is the lowest.

Usage

```

rain_reclus(cmeans.cluster)

## S3 method for class 'rain_reclus'
plot(x, ..., sample = "all", color.blind = FALSE)

```

Arguments

cmeans.cluster	an object of class cmeans_clus, inherited from cmeans_clus .
x	an object of class rain_reclus
...	Arguments to be passed to methods
sample	'all' to show all samples, or a numeric vector indicating the row number of samples in the sample table.
color.blind	logical. If TRUE colors optimized for colorblind readers are used.

Value

An object of class rain_reclus containing a sublist for each sample. Each sublist has the following components:

quality	quality threshold used in read_sample .
reference	reference ID.
centers	a data frame with the centroids coordinates.
data	a data frame with the fluorescence intensities and clusters name.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata",package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
samp <- read_sample(sample.table, system = "bio-rad",
                   file.location = fileLoc)

#Reference DBSCAN clustering
dbref <- reference_dbscan(ref, sample.table, save.template = FALSE)

#Predict position of clusters centroid from reference DBSCAN results
cent <- centers_data(samp, sample.table,dbref)

#Fuzzy c-means clustering
cmclus <- cmeans_clus(cent)

#Rain classification.
rainclus <- rain_reclus(cmclus)

plot(rainclus, sample = "all")

## End(Not run)
```

read_reference

Read reference files

Description

This function reads the results files of reference samples listed in the sample table. Fluoresce intensity and quality value (just for Thermo Fisher) are collected. If a [reference_dbscan](#) template file with the same input paramters (reference ID, eps, minPts) is available, fluorescence data, quality value and dbscan analysis results are retrived from the template file.

Usage

```
read_reference(
  sample.table,
  system = NULL,
  file.location = ".",
  reference.quality = 0.5,
  eps = NULL,
  minPts = NULL
)
```

Arguments

<code>sample.table</code>	object of class <code>sample_table</code> , inherited from read_sampleTable .
<code>system</code>	character. The name of digital PCR system used to generate the data. It must be either Thermo Fisher or Bio-Rad. Abbreviations are also accepted.
<code>file.location</code>	character. Full path name to reference and sample files location. The default corresponds to the working directory, (getwd). Tilde expansion (see (path.expand)) is performed.
<code>reference.quality</code>	numeric. Between 0 and 1. Quality threshold to subset the data. If different thresholds have to be applied to various reference samples, a vector of the same length of number of reference samples has to be provided. Used only when the system is Thermo Fisher.
<code>eps, minPts</code>	numeric. Input parameters for the DBSCAN algorithm. If they match the parameters of reference_dbscan template file, the data are retrieved from the template.

Value

An object of class `read_reference` containing a sublist for each reference. Each sublist has the following components:

<code>quality</code>	value of the <code>reference.quality</code> parameter.
<code>data</code>	a matrix with the fluorescence intensities and quality values.
<code>dbscan</code>	an object of class <code>dbscan_fast</code> , inherited from dbscan . This component is available only if a reference_dbscan template file is used to retrieve the data.

Examples

```
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
```

```

sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

```

read_sample	<i>Read sample files</i>
-------------	--------------------------

Description

This function reads the results files of samples listed in the sample table. Fluoresce intensity and quality value (just for Thermo Fisher) are collected.

Usage

```

read_sample(
  sample.table,
  system = NULL,
  file.location = ".",
  sample.quality = 0.5,
  partition.volume = NULL
)

```

Arguments

sample.table	object of class sample_table, inherited from read_sampleTable .
system	character. The name of digital PCR system used to generate the data. It must be either Thermo Fisher or Bio-Rad. Abbreviations are also accepted.
file.location	character. Full path name to reference and sample files location. The default corresponds to the working directory, (getwd). Tilde expansion (see (path.expand)) is performed.
sample.quality	numeric. Between 0 and 1. Quality threshold to subset data. If different thresholds have to be applied to various samples, a vector of the same length of number of samples has to be provided. Used only when the system is Thermo Fisher.
partition.volume	numeric. This parameters is taken into account when the parameter 'system' is set on Other. Indicate the partion volume in microliters spcific to the digital PCR system.

Value

An object of class read_sample containing a sublist for each sample. Each sublist has the following components:

quality	value of the sample.quality parameter.
data	a matrix with the fluorescence intensities and quality values.

Examples

```
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
samp <- read_sample(sample.table, system = "bio-rad",
                   file.location = fileLoc)
```

read_sampleTable	<i>Read sample table</i>
------------------	--------------------------

Description

This function reads a file containing the essential information about the samples and experimental settings. The file has to be filled out by the user and formatted as described in the vignette.

Usage

```
read_sampleTable(file, system = NULL, file.location = ".")
```

Arguments

file	character. The name or the path of csv file to be read. If it does not contain an absolute path, the file name is relative to the current working directory, (getwd).
system	character. The name of digital PCR system used to generate the data. It must be either Thermo Fisher or Bio-Rad. Abbreviations are also accepted.
file.location	character. Full path name to reference and sample files location. The default corresponds to the working directory, (getwd). Tilde expansion (see (path.expand)) is performed.

Value

An object of class `sample_table`.

Examples

```
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)
```

reference_dbscan	<i>Find the empty partitions and single target clusters in the reference sample</i>
------------------	---

Description

This function computes a DBSCAN analysis to identify single target clusters in the reference samples listed in the sample table. If a [reference_dbscan](#) template file with the same input parameters (reference ID, eps, minPts) is available, data are retrieved from the template file.

Usage

```
reference_dbscan(
  reference.subquality,
  sample.table,
  eps = 200,
  minPts = 50,
  save.template = FALSE
)

## S3 method for class 'reference_dbscan'
plot(x, ..., reference = "all")
```

Arguments

reference.subquality	an object of class read_reference, inherited from read_reference .
sample.table	object of class sample_table, inherited from read_sampleTable .
eps, minPts	numeric. Input parameters for the DBSCAN algorithm. If they match the parameters of reference_dbscan template file, the data are retrieved from the template.
save.template	logical. If TRUE a template of DBSCAN analysis of reference samples is saved. When system is Thermo Fisher, save.template can be also a character vector indicating the chipID.

x	an object of class reference_dbscan
...	Arguments to be passed to methods
reference	'all' to show all reference samples, or a character vector with chip ID (Thermo Fisher) or the file name (Bio-rad) of reference samples to be showed.

Value

An object of class reference_dbscan containing a sublist for each reference. Each sublist has the following components:

quality	quality threshold used in read_reference .
data	a matrix with the fluorescence intensities and quality values.
dbscan	an object of class dbscan_fast, inherited from dbscan .

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
samp <- read_sample(sample.table, system = "bio-rad",
                   file.location = fileLoc)

#Reference DBSCAN clustering
dbref <- reference_dbscan(ref, sample.table, save.template = FALSE)

plot(dbref, reference = "all")

## End(Not run)
```

replicates_quant	<i>Calculation of targets concentration, pooling the sample replicates</i>
------------------	--

Description

This function calculates the concentration of the targets, combining the results of the replicates of each sample.

Usage

```
replicates_quant(raw.results, sample.table)
```

Arguments

raw.results	an object of class target_quant, inherited from target_quant .
sample.table	object of class sample_table, inherited from read_sampleTable .

Value

An object of class replicates_quant containing a sublist for every sample. Each sublist has the following components:

quality	quality threshold used in read_sample .
reference	reference ID.
raw results	a data frame with the results of quantification.
replicates results	a data frame with the results of quantification of pooled replicates.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
```

```

samp <- read_sample(sample.table, system = "bio-rad",
                    file.location = fileLoc)

#Reference DBSCAN clustering
dbref <- reference_dbscan(ref, sample.table, save.template = FALSE)

#Predict position of clusters centroid from reference DBSCAN results
cent <- centers_data(samp, sample.table,dbref)

#Fuzzy c-means clustering
cmclus <- cmeans_clus(cent)

#Rain classification.
rainclus <- rain_reclus(cmclus)

#Quantification
quantcm <- target_quant(cmclus, sample.table)
quant <- target_quant(rainclus, sample.table)

#Replicates pooling
rep.quant <- replicates_quant(quant, sample.table)

## End(Not run)

```

report_dPCP

Export dPCP analysis results to a pdf report

Description

This function generates a pdf report of the dPCP analysis.

Usage

```
report_dPCP(data, filename, sample = "all", color.blind = FALSE)
```

Arguments

data	an object of class dPCP, inherited from dPCP .
filename	character. File name (no extension) for csv and pdf files to create on disk.
sample	'all' to show all samples, or a numeric vector indicating the row number of samples in the sample table.
color.blind	logical. If TRUE colors optimized for colorblind readers are used.

Value

A pdf file with the information and results of the dPCP analysis.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#dPCP analysis
results <- dPCP(sampleTable, system = "bio-rad", file.location = fileLoc,
               eps = 200, minPts = 50, save.template = FALSE,
               rain = TRUE)

report_dPCP(results, filename = "dPCRproject_1")

## End(Not run)
```

target_quant	<i>Calculation of targets concentration.</i>
--------------	--

Description

This function calculates the concentration of the targets according to the Poisson distribution.

Usage

```
target_quant(data.cluster, sample.table)
```

Arguments

data.cluster	an object of class rain_reclus or cmeans_clus.
sample.table	object of class sample_table, inherited from read_sampleTable .

Value

An object of class target_quant containing a sublist for each sample. Each sublist has the following components:

quality	quality threshold used in read_sample .
reference	reference ID.
raw results	a data frame with the results of the quantification.

Examples

```
## Not run:
library(dPCP)

#Find path of sample table and location of reference and input files
sampleTable <- system.file("extdata", "Template_sampleTable.csv",
                           package = "dPCP")

fileLoc <- system.file("extdata", package = "dPCP")

#Read sample table file
sample.table <- read_sampleTable(sampleTable, system = "bio-rad",
                                file.location = fileLoc)

#Read reference files
ref <- read_reference(sample.table, system = "bio-rad",
                     file.location = fileLoc)

#Read samples files
samp <- read_sample(sample.table, system = "bio-rad",
                   file.location = fileLoc)

#Reference DBSCAN clustering
dbref <- reference_dbscan(ref, sample.table, save.template = FALSE)

#Predict position of clusters centroid from reference DBSCAN results
cent <- centers_data(samp, sample.table,dbref)

#Fuzzy c-means clustering
cmclus <- cmeans_clus(cent)

#Rain classification.
rainclus <- rain_reclus(cmclus)

#Quantification
quantcm <- target_quant(cmclus, sample.table)
quant <- target_quant(rainclus, sample.table)

## End(Not run)
```

Index

centers_data, [2](#), [4](#)
cmeans, [4](#)
cmeans_clus, [3](#), [11](#)

dbscan, [6](#), [7](#), [13](#), [17](#)
dbscan_combination, [5](#)
dPCP, [6](#), [10](#), [19](#)

export_csv, [9](#)

getwd, [5](#), [7](#), [13–15](#)

manual_correction, [10](#)

path.expand, [5](#), [7](#), [13–15](#)
plot.centers_data (centers_data), [2](#)
plot.cmeans_clus (cmeans_clus), [3](#)
plot.dPCP (dPCP), [6](#)
plot.rain_reclus (rain_reclus), [11](#)
plot.reference_dbscan
 (reference_dbscan), [16](#)

rain_reclus, [11](#)
read_reference, [12](#), [16](#), [17](#)
read_sample, [2–4](#), [11](#), [14](#), [18](#), [20](#)
read_sampleTable, [2](#), [13](#), [14](#), [15](#), [16](#), [18](#), [20](#)
reference_dbscan, [2](#), [12](#), [13](#), [16](#), [16](#)
replicates_quant, [18](#)
report_dPCP, [19](#)

target_quant, [18](#), [20](#)